

-----  
\* Introduction  
-----

"AutoInt1.dll" plug-in (32-bit) implements COM automation programming interface (ProgId "EchoWave2.CmdInt1") that allows to access "Echo Wave II" ultrasound images from other applications and adjust "Echo Wave II" ultrasound scanning parameters.

In order to use this interface, in "Echo Wave II" file "...\\Config\\Plugins\\list.txt" must be present the following line:

AutoInt1.dll

It means that during software startup will be loaded plug-in dll.

In order to use automation interface, "Echo Wave II" software must be running.

"Echo Wave II" software is run using Administrator rights. This means that another application that will be connecting to "Echo Wave II" also must be run As Administrator.

"AutoInt1Client.dll" is an "Echo Wave II" automation interface client .Net assembly that can be used from 32/64-bit applications that support .Net.

"AutoInt1Client.cs" - C# source code of "AutoInt1Client.dll" that can be used as a reference for implementing automation interface clients in other programming languages.

"AutoInt1Client\_test.m" - MATLAB (R2010b or more recent version) script that demonstrates how to use .Net client dll and access "Echo Wave II" ultrasound images. In order to test the script, please remove/add comments (%) from/to desired script lines.

Programming interface allows to do the following things: freeze/run or record/stop ultrasound (when is connected ultrasound scanner), open TVD or other file that is supported by "Echo Wave II", get the number of cine frames, select desired frame, get frame time in milliseconds (the time of first frame is 0.0), get frame in grayscale or RGB format, change ultrasound scanning mode, set ultrasound scanning parameters.

IMPORTANT. If some command opens "Echo Wave II" modal dialog (window) like save, open, print, patient window, options window or any message box, this means that called software will be blocked until the user will close that dialog in "Echo Wave II" software. The same applies to all other commands. If you want to send commands to "Echo Wave II" and do not block user interface of your software, then connection to "Echo Wave II" must be done and commands must be sent from a single separate thread.

-----  
\* AutoInt1Client.dll functions  
-----

Namespace: AutoInt1Client

Class: CmdInt1

In order to use below listed functions, must be created CmdInt1 object. Since these functions call appropriate "Echo Wave II" functions, their operation depends on what is enabled in "Echo Wave II" user interface.

-----  
int ConnectToRunningProgram()  
-----

Connects to "Echo Wave II" automation interface with ProgId "EchoWave2.CmdInt1". If connection fails, this function returns "-1". In other case it returns "0". Connection is possible only if "Echo Wave II" is running and both connected applications were run "As Administrator".

-----  
int OpenFile(string fn)  
-----

Open file by "Echo Wave II".  
fn - full path to file, for example "C:\\Echo Images\\test1.tvd"  
In case of unexpected errors this function returns "-1".

-----  
int GetFramesCount()  
-----

Get the number of frames in opened file or scanned cine.  
In case of unexpected errors this function returns "-1".

-----  
int GoToFrame1n(int frm\_idx\_1n, bool load\_frame\_data)  
-----

Go to frame with index frm\_idx\_1n. Indexing starts from "1".  
load\_frame\_data must be true if we will want to get image (bytes) of selected frame using functions GetLoadedFrame\*.  
Before using this function, must be stopped cine playback, stopped recording, stopped ultrasound scanning.  
If, although not recommended, this function will be used without stopping ultrasound playback/recording/scanning in order to load current frame, then parameter frm\_idx\_1n must be "-1".  
In case of unexpected errors this function returns "-1".

-----  
int GetCurrentFrameIdx1n()  
-----

Get index of current (visible in "Echo Wave II") frame.  
In case of unexpected errors this function returns "-1".

-----  
double GetCurrentFrameTime()  
-----

Get time (in milliseconds) of current frame.  
The time of first frame is 0.0.  
In case of unexpected errors this function returns "-1".

-----  
int GetLoadedFrameWidth()  
-----

Get loaded frame width in pixels.  
In case of unexpected errors this function returns "-1".

-----  
int GetLoadedFrameHeight()  
-----

Get loaded frame height in pixels.  
In case of unexpected errors this function returns "-1".

-----  
int[, ,] GetLoadedFrameRGB()  
-----  
Get three-dimensional array of loaded frame pixel data.  
Returned array contains red (R), green (G), blue(B) planes of image pixels.  
R, G, B values are in the interval [0;255].  
In case of unexpected errors this function returns null.

-----  
int[,] GetLoadedFrameGray()  
-----  
Get two-dimensional array of loaded frame grayscale pixel data.  
Returned grayscale  $((R+G+B)/3)$  values are in the interval [0;255].  
In case of unexpected errors this function returns null.

-----  
int[,] GetLoadedFrameRed()  
-----  
Get two-dimensional array of Red image (24-bit RGB) plane.  
Returned values are in the interval [0;255].  
In case of unexpected errors this function returns null.

-----  
int[,] GetLoadedFrameGreen()  
-----  
Get two-dimensional array of Green image (24-bit RGB) plane.  
Returned values are in the interval [0;255].  
In case of unexpected errors this function returns null.

-----  
int[,] GetLoadedFrameBlue()  
-----  
Get two-dimensional array of Blue image (24-bit RGB) plane.  
Returned values are in the interval [0;255].  
In case of unexpected errors this function returns null.

-----  
int GetUltrasoundX1(int img\_id)  
-----  
Get X1 coordinate in pixels of ultrasound image region with passed img\_id.  
Here img\_id: 1-B; 2-B2; 3-B3; 4-B4; 7-M; 8-PW; 9-CW.  
In case of unexpected errors this function returns "-1".

-----  
int GetUltrasoundX2(int img\_id)  
-----  
Get X2 coordinate in pixels of ultrasound image region with passed img\_id.  
Here img\_id: 1-B; 2-B2; 3-B3; 4-B4; 7-M; 8-PW; 9-CW.  
In case of unexpected errors this function returns "-1".

-----  
int GetUltrasoundY1(int img\_id)  
-----  
Get Y1 coordinate in pixels of ultrasound image region with passed img\_id.  
Here img\_id: 1-B; 2-B2; 3-B3; 4-B4; 7-M; 8-PW; 9-CW.  
In case of unexpected errors this function returns "-1".

-----  
int GetUltrasoundY2(int img\_id)  
-----

Get Y2 coordinate in pixels of ultrasound image region with passed img\_id.  
Here img\_id: 1-B; 2-B2; 3-B3; 4-B4; 7-M; 8-PW; 9-CW.  
In case of unexpected errors this function returns "-1".

-----  
double GetUltrasoundPhysicalDeltaX(int img\_id)  
-----

Get horizontal resolution of ultrasound image region with passed img\_id.  
Here img\_id: 1-B; 2-B2; 3-B3; 4-B4; 7-M; 8-PW; 9-CW.  
Returned value is in cm/pt (centimeters per pixel) or s/pt (seconds per pixel)  
depending on image type.  
In case of unexpected errors this function returns "0".

-----  
double GetUltrasoundPhysicalDeltaY(int img\_id)  
-----

Get vertical resolution of ultrasound image region with passed img\_id.  
Here img\_id: 1-B; 2-B2; 3-B3; 4-B4; 7-M; 8-PW; 9-CW.  
Returned value is in cm/pt or cm/s/pt depending on image type.  
In case of unexpected errors this function returns "0".

-----  
int FreezeRun()  
-----

Freeze/Run ultrasound scanning.  
In case of unexpected errors this function returns "-1".

-----  
int IsRunState()  
-----

Returns "1" if ultrasound scanning is in progress.  
In case of unexpected errors this function returns "-1".

-----  
int PlayPause()  
-----

Play/Pause scanned or opened cine.  
This function works when is enabled "Echo Wave II" Play/Pause button.  
In case of unexpected errors this function returns "-1".

-----  
int IsPlayState()  
-----

Returns "1" if cine playback is in progress.  
In case of unexpected errors this function returns "-1".

-----  
int RecordStop()  
-----

Record/Stop ultrasound cine (when ultrasound scanner is connected).  
When recording is started, software clears cine buffer.  
When recording is stopped, software freezes ultrasound scanning.

This function works only when is enabled "Echo Wave II" Record/Stop button. That is, at first ultrasound scanning must be run using FreezeRun(), then found anatomy of interest, then called RecordStop() in order to clear cine and do recording and then called another RecordStop() in order to stop recording. In case of unexpected errors this function returns "-1".

```
-----  
int IsRecordingState()  
-----
```

Returns "1" if cine recording is in progress.  
In case of unexpected errors this function returns "-1".

```
-----  
int WmCopyDataCmd(string cmd_str)  
-----
```

Invoke WM\_COPYDATA command cmd\_str. For the list of supported WM\_COPYDATA commands please check "Echo Wave II" user manual section "Software control using command line".  
In case of unexpected errors this function returns "-1".

```
-----  
int ParamSet(int param_id, int val)  
-----
```

Shift parameter value by specified index val or set parameter value val depending on passed param\_id. The list of defines of param\_id is presented below (#define id...).  
If we have B Gain values 10%, 20%, 30%, 40%, 50%, current value is 40%, then param\_id=id\_b\_gain\_shift and val=-1 will set Gain value to 30%.  
If some command does not need parameters, then should be passed val=0.  
In case of unexpected errors this function returns "-1".

For cases when it is not enough to pass single int value or when defined command returns some value, are used special functions that will be mentioned near each define.  
These functions can be used only in cases where for concrete param\_id is documented that must be used such functions. In other cases you will get either exceptions or crashes because of attempt to pass/return incorrect data type(s).

```
int ParamSet2(int param_id, int val, int val2)  
int ParamSetString(int param_id, string val)  
object ParamGet(int param_id)  
int ParamGetInt(int param_id)  
int ParamGetInt2(int param_id, int val)  
bool ParamGetBool(int param_id)  
bool ParamGetBool2(int param_id, int val)  
double ParamGetDouble(int param_id)  
float ParamGetFloat(int param_id)  
string ParamGetString(int param_id)  
string ParamGetString2(int param_id, int val)  
string ParamGetString3(int param_id, int val, val2)
```

Below is the list of supported param\_id defines and their description:

```
#define id_get_current_beamformer_code 915 // GET current beamformer code  
as int. ParamGetInt(id_get_current_beamformer_code)
```

```

#define id_get_current_beamformer_name 916 // GET current beamformer name
as string. ParamGetString(id_get_current_beamformer_name)
#define id_get_current_probe_code 917 // GET current probe code as
int. ParamGetInt(id_get_current_probe_code)
#define id_get_current_probe_name 918 // GET current probe name as
string. ParamGetString(id_get_current_probe_name)
#define id_probe 104 // Probe button; val = 0;
#define id_is_no_probe 223 // Get if the probe is not detected.
ParamGetBool(id_is_no_probe)
#define id_is_unsupported_probe 224 // Get if is connected unsupported
probe. ParamGetBool(id_is_unsupported_probe)
#define id_is_probe_active 103 // Get if the probe is active (the probe
becomes inactive when is opened tpd/tvd file). ParamGetBool(id_is_probe_active)
#define id_get_probes_count 225 // Get the number of connected probes.
ParamGetInt(id_get_probes_count)
#define id_freeze_run 100 // Freeze/Run command; val = 0;
#define id_is_freeze 101 // GET if ultrasound scanning is frozen.
ParamGetBool(id_is_freeze)
#define id_cine_play_pause 119 // Cine play/pause; val = 0;
#define id_cine_pause 141 // Cine pause button; val = 0;
#define id_cine_record_stop 142 // Cine record/stop button; The image
must be unfrozen. First invocation clears cine buffers, second invocation does
freeze. val = 0;
#define id_cine_frame_shift 120 // Cine previous or next frame
#define id_is_cine_playing 122 // GET from software if cine is
playing. ParamGetBool(id_is_cine_playing)
#define id_is_cine_record_in_progress 143 // Get if cine is in recording
mode (not simple scanning, but recording).
ParamGetBool(id_is_cine_record_in_progress)
#define id_is_usg_file_opened 191 // Returns if we have opened
tpd/tvd file. ParamGetBool(id_is_usg_file_opened)
#define id_auto_adjust 116 // Automatic image adjustment (Auto
button); val = 0;

// Scanning modes
#define id_button_b 174 // B button's command; val = 0;
#define id_button_dual_quad 168 // Dual/Quad button's last used
command; val = 0;
#define id_button_dual 898 // B+B (Dual) button; val = 0;
#define id_button_quad 899 // 4B (Quad) button; val = 0;
#define id_button_m 900 // B+M button; val = 0;
#define id_button_cd 901 // CD button's last used
command; val = 0;
#define id_button_pw 625 // PW button (mode); val = 0;
#define id_button_cw 662 // CW button (mode); val = 0;
#define id_button_pw_cw 897 // PW/CW button's last used
command; val = 0;
#define id_button_update 902 // "Update" (F4) button; val =
0;
#define id_use_duplex_triplex_shift 871 // Option "Use Duplex/Triplex
modes"; -1 - off; 1 - on;
#define id_use_duplex_triplex 895 // Option "Use Duplex/Triplex
modes"; toggle; val = 0;
#define id_is_triplex_supported 227 // GET if is supported Triplex
(simultaneous running B+CFM+PW). ParamGetBool(id_is_triplex_supported)
#define id_get_use_duplex_triplex 228 // GET the value of option "Use

```

```

Duplex/Triplex modes". ParamGetBool(id_get_use_duplex_triplex)
#define id_is_color_m_supported      247      // GET if is supported Color M
mode. ParamGetBool(id_is_color_m_supported)

// Set scanning state. This command can be used for cases that are not available
as
// commands of scanning mode buttons.
// As parameter is passed scanning state id (e.g., id_state_b_b1r).
// In many cases the sequence of states is important, so before programming it
is recommended
// to check in "Echo Wave II" how "Echo Wave II" buttons change scanning states.
// If direct change from one state to another is unsupported by "Echo Wave II",
new state
// will be set by at first entering id_state_b_b1r.
// In state identifier names "r" means running, "f" means frozen.
#define id_scanning_state_set      201      // SET scanning state.
ParamSet(id_scanning_state_set, id_state_b_b1r)
#define id_scanning_state_get      200      // GET current scanning state.
ParamGetInt(id_kb_cmd_scanning_state_get)

// B mode states
#define id_state_b_b1r              1
#define id_state_b_b1f              2

// B+M, M states
#define id_state_bm_b1mr            3
#define id_state_bm_b1mf            4
#define id_state_bm_m1r            5
#define id_state_bm_m1f            6
#define id_state_bm_b1mr_m1r       7
#define id_state_bm_b1mf_m1f       8

// Dual states (B+B states)
#define id_state_bb_b1r_b2f         11
#define id_state_bb_b1f_b2r         13
#define id_state_bb_b1f             14
#define id_state_bb_b2f             15
#define id_state_bb_b1f_b2f_sel_b1  22
#define id_state_bb_b1f_b2f_sel_b2  23
#define id_state_bb_b1f_b2f         26
#define id_state_bb_b1r_b2r         27

// Color Doppler states (CFM, PDI, DPDI)
#define id_state_cfm_cfm1r          16
#define id_state_cfm_cfm1f          17
#define id_state_pdi_pdi1r          18
#define id_state_pdi_pdi1f          19
#define id_state_dpdi_dpdi1r        20
#define id_state_dpdi_dpdi1f        21

// Quad states (4B states)
#define id_state_4b_b1r_b2f_b3f_b4f 30
#define id_state_4b_b1f_b2r_b3f_b4f 31
#define id_state_4b_b1f_b2f_b3r_b4f 32
#define id_state_4b_b1f_b2f_b3f_b4r 33
#define id_state_4b_b1f_b2f_b3f_b4f_sel_b1 34

```

```

#define id_state_4b_b1f_b2f_b3f_b4f_sel_b2      35
#define id_state_4b_b1f_b2f_b3f_b4f_sel_b3      36
#define id_state_4b_b1f_b2f_b3f_b4f_sel_b4      37

// Duplex states (B+PW)
#define id_state_bpw_b1pwr                        40 // B+(PW line) running
#define id_state_bpw_b1pwf                        41
#define id_state_bpw_pw1r                        42 // PW running
#define id_state_bpw_pw1f                        43
#define id_state_bpw_b1pwr_pw1r                 44 // B+PW runing
#define id_state_bpw_b1pwf_pw1f                 45

// Triplex states (CFM+PW)
#define id_state_cfm1pwr                         46
#define id_state_cfm1pwf                         47
#define id_state_cfm1pwr_pw1r                   48
#define id_state_cfm1pwr_pw1f                   49
#define id_state_cfm1pwr_pw1r_pw1r              50
#define id_state_cfm1pwr_pw1f_pw1r              51

// Triplex states (PDI+PW)
#define id_state_pdipw_pdi1pwr                   52
#define id_state_pdipw_pdi1pwf                   53
#define id_state_pdipw_pw1r                       54
#define id_state_pdipw_pw1f                       55
#define id_state_pdipw_pdi1pwr_pw1r              56
#define id_state_pdipw_pdi1pwr_pw1f              57

// Triplex states (DPDI+PW)
#define id_state_dpdi1pwr                         58
#define id_state_dpdi1pwf                         59
#define id_state_dpdi1pwr_pw1r                   60
#define id_state_dpdi1pwr_pw1f                   61
#define id_state_dpdi1pwr_pw1r_pw1r              62
#define id_state_dpdi1pwr_pw1f_pw1r              63

// Duplex Update states (B+PW)
#define id_state_bpw_b1pwr_pw1f                 64 // B running + PW frozen
#define id_state_bpw_b1pwr_pw1f_pw1r            65 // B frozen + PW running
#define id_state_bpw_b1pwr_pw1f_sel_b           66 // B frozen + PW frozen and selected
is B
#define id_state_bpw_b1pwr_pw1f_sel_pw         67 // B frozen + PW frozen and selected
is PW

// "Triplex" Update states (CFM+PW)
#define id_state_cfm1pwr_pw1f                   68 // CFM running + PW frozen
#define id_state_cfm1pwr_pw1f_pw1r              69 // CFM frozen + PW running
#define id_state_cfm1pwr_pw1f_sel_cfm           70 // CFM frozen + PW frozen
and selected is CFM
#define id_state_cfm1pwr_pw1f_sel_pw           71 // CFM frozen + PW frozen
and selected is PW

// "Triplex" Update states (PDI+PW)
#define id_state_pdipw_pdi1pwr_pw1f             72 // PDI running + PW frozen
#define id_state_pdipw_pdi1pwr_pw1f_pw1r        73 // PDI frozen + PW running
#define id_state_pdipw_pdi1pwr_pw1f_sel_pdi     74 // PDI frozen + PW frozen

```



```

and selected is PDI
#define id_state_pdipw_pdi1pwf_pw1f_sel_pw 75 // PDI frozen + PW frozen
and selected is PW

// "Triplex" Update states (DPDI+PW)
#define id_state_dpdiw_dpdi1pwr_pw1f 76 // DPDI running + PW frozen
#define id_state_dpdiw_dpdi1pwf_pw1r 77 // DPDI frozen + PW running
#define id_state_dpdiw_dpdi1pwf_pw1f_sel_dpdi 78 // DPDI frozen + PW frozen
and selected is DPDI
#define id_state_dpdiw_dpdi1pwf_pw1f_sel_pw 79 // DPDI frozen + PW frozen
and selected is PW

// B+CW (Continuous Wave Doppler)
#define id_state_bcw_b1cwr 97 // B+(CW line) running
#define id_state_bcw_b1cwf 98
#define id_state_bcw_cw1r 99 // CW running
#define id_state_bcw_cw1f 100
#define id_state_bcw_b1cwr_cw1f 103 // B running + CW frozen
#define id_state_bcw_b1cwf_cw1r 104 // B frozen + CW running
#define id_state_bcw_b1cwf_cw1f_sel_b 105 // B frozen + CW frozen and
selected is B
#define id_state_bcw_b1cwf_cw1f_sel_cw 106 // B frozen + CW frozen and
selected is CW

// CFM+CW
#define id_state_cfmw_cfm1cwr 107
#define id_state_cfmw_cfm1cwf 108
#define id_state_cfmw_cw1r 109
#define id_state_cfmw_cw1f 110
#define id_state_cfmw_cfm1cwr_cw1f 113 // CFM running + CW frozen
#define id_state_cfmw_cfm1cwf_cw1r 114 // CFM frozen + CW running
#define id_state_cfmw_cfm1cwf_cw1f_sel_cfm 115 // CFM frozen + CW frozen
and selected is CFM
#define id_state_cfmw_cfm1cwf_cw1f_sel_cw 116 // CFM frozen + CW frozen
and selected is CW

// PDI+CW
#define id_state_pdicw_pdi1cwr 117
#define id_state_pdicw_pdi1cwf 118
#define id_state_pdicw_cw1r 119
#define id_state_pdicw_cw1f 120
#define id_state_pdicw_pdi1cwr_cw1f 123 // PDI running + CW frozen
#define id_state_pdicw_pdi1cwf_cw1r 124 // PDI frozen + CW running
#define id_state_pdicw_pdi1cwf_cw1f_sel_pdi 125 // PDI frozen + CW frozen
and selected is PDI
#define id_state_pdicw_pdi1cwf_cw1f_sel_cw 126 // PDI frozen + CW frozen
and selected is CW

// DPDI+CW
#define id_state_dpdiw_dpdi1cwr 127
#define id_state_dpdiw_dpdi1cwf 128
#define id_state_dpdiw_cw1r 129
#define id_state_dpdiw_cw1f 130
#define id_state_dpdiw_dpdi1cwr_cw1f 133 // DPDI running + CW frozen
#define id_state_dpdiw_dpdi1cwf_cw1r 134 // DPDI frozen + CW running
#define id_state_dpdiw_dpdi1cwf_cw1f_sel_dpdi 135 // DPDI frozen + CW frozen

```

```

and selected is DPDI
#define id_state_dpdicw_dpdi1cwf_cw1f_sel_cw    136 // DPDI frozen + CW frozen
and selected is CW

// B+CM, CM states (Color M; M Color Flow Mode)
#define id_state_bcm_cfm1mr                    137 // M line on CFM image
#define id_state_bcm_cfm1mf                    138
#define id_state_bcm_cm1r                      139 // pure Color M
#define id_state_bcm_cm1f                      140
#define id_state_bcm_b1mr_cm1r                 141 // B + Color M
#define id_state_bcm_b1mf_cm1f                 142

// B mode controls
#define id_b_frequency_shift                    300 // B Frequency
#define id_b_frequency_shift_loop              320 // B Frequency; go to
first when is reached last value
#define id_b_change_thi                        178 // B Switch THI on/off;
val = 0;
#define id_b_is_thi_frequency                  177 // GET if current
frequency is THI. ParamGetBool(id_b_is_thi_frequency)
#define id_b_focus_shift                       302 // B Focus depth
#define id_b_focus_shift_loop                 321 // B Focus depth; go to
first when is reached last value
#define id_b_focuses_number_shift             334 // B Focuses Number
#define id_b_focus_set_shift                   335 // B Focus Set
#define id_b_dynamic_focus                     304 // B Dynamic Focus
on/off; val = 0;
#define id_b_dynamic_focus_shift              323 // B Dynamic Focus; pass
negative or positive value in order to turn on or off
#define id_b_is_dynamic_focus                  171 // GET if is turned on
Dynamic Focus. ParamGetBool(id_b_is_dynamic_focus)
#define id_b_depth_shift                       305 // B Depth
#define id_b_depth_shift_asynch               306 // B Depth asynchronous
shift
#define id_b_power_shift                       307 // B Power
#define id_b_gain_shift                        309 // B Gain
#define id_b_dynamic_range_shift              311 // B Dynamic Range
#define id_b_palette_shift                     338 // B Palette (B Color
Map)
#define id_b_palette_gamma_shift              313 // B Palette Gamma
#define id_b_palette_brightness_shift         315 // B Palette Brightness
#define id_b_palette_contrast_shift           317 // B Palette Contrast
#define id_b_palette_negative                  319 // B Palette Negative;
val = 0;
#define id_b_dynamic_range_shift_loop          322 // B Dynamic Range; go
to first when is reached last value
#define id_b_view_area_shift                   324 // B View Area
#define id_b_view_area_shift_loop             325 // B View Area; go to
first when is reached last value
#define id_b_frame_averaging_shift            326 // B Frame Averaging
#define id_b_rejection_shift                   327 // B Rejection
#define id_b_image_enhancement_enabled_shift  328 // B Image Enhancement
Enabled; negative value - turn off; positive - turn on
#define id_b_image_enhancement_method_shift   329 // B Image Enhancement
Method
#define id_b_image_enhancement_method_shift2  336 // B Image Enhancement;

```

```

Off, 1, 2, ...
#define id_b_speckle_reduction_enabled_shift    330    // B Speckle Reduction
Enabled; negative value - turn off; positive - turn on
#define id_b_speckle_reduction_level_shift      331    // B Speckle Reduction
Level
#define id_b_speckle_reduction_level_shift2     337    // B Speckle Reduction:
Off, 1, 2, ...
#define id_b_lines_density_shift                332    // B Lines Density
#define id_b_change_scan_direction             106    // B Change scan
direction; val = 0;
#define id_b_is_scan_direction_changed          133    // GET bool scan
direction value of current image. ParamGetBool(id_b_is_scan_direction_changed)
#define id_b_flip_up_down                       105    // B Flip image up/down;
val = 0;
#define id_b_rotate_shift                       333    // B Rotate
#define id_b_rotate_get                         132    // GET int rotation
value of current image. ParamGetInt(id_b_rotate_get)
#define id_b_steering_trapezoid_angle_shift     339    // B Angle
(Steering/Trapezoid)
#define id_b_steering_trapezoid_angle_shift_loop 341    // B Angle
(Steering/Trapezoid); go to first when is reached last value
#define id_b_scan_type_shift                    340    // B Scan Type:
Standard, Trapezoid, Compound
#define id_b_scan_type_shift_loop              344    // B Scan Type:
Standard, Trapezoid, Compound; go to first when is reached last value
#define id_b_compound_on_off                    342    // B Compound on/off;
val = 0;
#define id_zoom_controls                        110    // Show Zoom ControlBar;
val = 0;
#define id_b_zoom_default                       111    // B Set zoom 1:1
(default); val = 0;
#define id_b_zoom_factor_shift                  112    // B Zoom
#define id_b_zoom_factor_shift_asynch          113    // B Zoom asynchronous
#define id_b_tgc_shift                          114    // B TGC; As parameter
is passed the number (index) of TGC control [0;4] and shift of value.
ParamSet2(id_b_tgc_shift, 0, 1);
#define id_b_tgc_get_min                        179    // GET minimal possible
TGC value. ParamGetInt(id_b_tgc_get_min)
#define id_b_tgc_get_max                        180    // GET maximal possible
TGC value. ParamGetInt(id_b_tgc_get_max)
#define id_b_tgc_set                            181    // Set TGC value. As
parameter is passed the number (index) of TGC control [0;4] and new value.
ParamSet2(id_b_tgc_set, 0, 10);

// M mode controls
#define id_m_line_position_shift                400    // M Line Position
#define id_m_line_position_shift_asynch         401    // M Line Position
asynchronous shift
#define id_m_sweep_speed_shift                  402    // M Sweep Speed
#define id_m_palette_gamma_shift                404    // M Palette Gamma
#define id_m_palette_brightness_shift           406    // M Palette Brightness
#define id_m_palette_contrast_shift             408    // M Palette Contrast
#define id_m_palette_negative                   410    // M Palette Negative;
val = 0;
#define id_m_palette_negative_shift             411    // M Palette Negative;
-1 - off, +1 - on

```

```

#define id_m_rejection_shift          412      // M Rejection
#define id_m_zoom_default            413      // M Zoom 1:1 (default);
val = 0;
#define id_m_zoom_factor_shift       414      // M Zoom factor
#define id_m_zoom_position_shift     415      // M Zoom vertical shift

// Color Doppler (CD) mode controls (CFM, PDI, DPDI)
#define id_cd_frequency_shift        500      // CD Frequency
#define id_cd_frequency_shift_loop   501      // CD Frequency; go to
first when is reached last value
#define id_cd_steering_angle_shift   502      // CD Steering Angle (CD
Angle)
#define id_cd_steering_angle_shift_loop 536    // CD Steering Angle (CD
Angle); go to first when is reached last value
#define id_cd_prf_shift              504      // CD PRF
#define id_cd_power_shift            506      // CD Power
#define id_cd_gain_shift             508      // CD Gain (changes CFM
or PDI/DPDI Gain depending on current scanning mode)
#define id_cd_window_horiz_position_shift 510   // CD Window Horizontal
Position
#define id_cd_window_horiz_position_shift_asynch 538    // CD Window
Horizontal Position asynchronous shift
#define id_cd_window_vert_position_shift 512    // CD Window Vertical
Position
#define id_cd_window_vert_position_shift_asynch 539    // CD Window
Vertical Position asynchronous shift
#define id_cd_window_size_shift      537      // CD Window Size
#define id_cd_window_size_shift_loop 513      // CD Window Size; go to
first when is reached last value
#define id_cd_window_horiz_size_shift 514      // CD Window Horizontal
Size
#define id_cd_window_vert_size_shift 516      // CD Window Vertical
Size
#define id_cd_window_horiz_vert_size_shift 540     // CD Window
Horizontal and Vertical Size
#define id_cd_window_horiz_vert_size_shift_asynch 541    // CD Window
Horizontal and Vertical Size asynchronous shift
#define id_cd_color_map_shift        518      // CD Color Map (CFM,
PDI, DPDI)
#define id_cd_palette_invert_shift   531      // CD Palette Invert
(CFM or DPDI depending on current scanning mode); -1 - invert on; 1 - invert
off;
#define id_cd_palette_invert        535      // CD Palette Invert
(CFM or DPDI depending on current scanning mode); Toggle on/off; val = 0;
#define id_cd_color_priority_shift   520      // CD B/Color Priority
#define id_cd_color_threshold_shift  522      // CD Color Threshold
#define id_cd_lines_density_shift    524      // CD Lines Density
#define id_cd_wall_filter_shift      526      // CD Wall Filter
#define id_cd_baseline_shift         527      // CD Baseline
#define id_cd_color_averaging_shift  529      // CD Color Averaging
#define id_cd_dynamic_range_shift    530      // PDI/DPDI Dynamic
Range
#define id_cd_pulses_number_shift    532      // CD Pulses Number (CD
Pulse Length)
#define id_cd_packet_size_shift      533      // CD Packet Size (CFM
or PDI/DPDI depending on current scanning mode)

```

```

#define id_cd_spatial_filtering_shift2      534    // CD Spatial Filtering:
Off, 1, 2, 3, 4
#define id_cd_scale_shift                   542    // PDI/DPDI Scale
#define id_cd_transparency_shift2          543    // CD Transparency: Off,
0.1, 0.2, ...

// Pulsed Wave (PW) and Continuous Wave (CW) Doppler mode controls
#define id_pw_correction_angle_shift        600    // PW Correction Angle
#define id_pw_steering_angle_shift          602    // PW Steering Angle
#define id_pw_steering_angle_shift_loop    641    // PW Steering Angle go to first
when is reached last value
#define id_pw_prf_shift                     604    // PW PRF
#define id_pw_color_map_shift              605    // PW Color Map
#define id_pw_palette_invert_shift         603    // PW Palette Invert; -1 - off;
1 - on;
#define id_pw_palette_invert               681    // PW Palette Invert; toggle;
val = 0;
#define id_pw_power_shift                   606    // PW Power
#define id_pw_dynamic_range_shift          607    // PW Dynamic Range
#define id_pw_gain_shift                   608    // PW Gain
#define id_pw_sample_volume_horiz_position_shift 610    // PW Sample
Volume Horizontal Position
#define id_pw_sample_volume_horiz_position_shift_asynch2 642    // PW
Sample Volume Horizontal Position
#define id_pw_sample_volume_vert_position_shift 612    // PW Sample
Volume Vertical Position
#define id_pw_sample_volume_vert_position_shift_asynch2 643    // PW
Sample Volume Vertical Position
#define id_pw_sample_volume_size_shift      614    // PW Sample
Volume Size
#define id_pw_sweep_speed_shift             616    // PW Sweep Speed
#define id_pw_baseline_shift               618    // PW Baseline
#define id_kb_cmd_pw_baseline_shift2       682    // PW Baseline shift taking into
account PW Invert
#define id_pw_invert                       620    // PW Invert checkbox; toggle;
val = 0;
#define id_pw_is_pw_invert                 134    // GET bool PW invert value of
current image. ParamGetBool(id_pw_is_pw_invert)
#define id_pw_sound_volume_shift           621    // PW Sound Volume
#define id_pw_wall_filter_shift            623    // PW Wall Filter
#define id_pw_invert_shift                 626    // PW Invert; -1 - off, +1 - on
#define id_pw_scale_shift                  627    // PW Scale
#define id_pw_frequency_shift              628    // PW Frequency
#define id_pw_frequency_shift_loop         629    // PW Frequency; go to first
when is reached last value
#define id_pw_smooth_shift                 639    // PW Smoothing
#define id_pw_spectral_averaging           640    // PW Spectral Averaging;
toggle; val = 0;
#define id_pw_auto_trace                   680    // Turn PW/CW auto trace on/off;
toggle; val = 0;

#define id_are_controls_enabled_b          202    // Return true if B controls are
enabled. ParamGetBool(id_are_controls_enabled_b)
#define id_are_controls_enabled_m          203    // Return true if M controls are
enabled. ParamGetBool(id_are_controls_enabled_m)
#define id_are_controls_enabled_pw         204    // Return true if PW controls

```

```

are enabled. ParamGetBool(id_are_controls_enabled_pw)
#define id_are_controls_enabled_cfm      205      // Return true if CFM controls
are enabled. ParamGetBool(id_are_controls_enabled_cfm)
#define id_are_controls_enabled_pdi      206      // Return true if PDI controls
are enabled. ParamGetBool(id_are_controls_enabled_pdi)
#define id_are_controls_enabled_dpdi     207      // Return true if DPDI controls
are enabled. ParamGetBool(id_are_controls_enabled_dpdi)

#define id_is_control_enabled            208      // Return true if control with
passed integer identifier is enabled. ParamGetBool2(id_is_control_enabled,
b_focus_shift)
// Supported id_is_control_enabled identifiers: m_line_position_shift,
cd_window_position_shift, pw_sample_volume_position_shift,
// b_focus_shift, zoom_factor_shift, pw_sweep_speed_shift, m_sweep_speed_shift,
b_steering_trapezoid_angle_shift
// cd_steering_angle_shift, pw_steering_angle_shift, cd_prf_shift,
cd_window_size_shift,
// pw_prf_shift, pw_scale_shift, use_duplex_triplex_shift, cd_gain_shift,
pw_gain_shift,
// pw_baseline_shift, thumbnails_show, cd_baseline_shift,
pw_correction_angle_shift, pw_sample_volume_position_shift.

#define id_get_exam_id                  236      // GET identifier of current
exam (from "Patient" window). ParamGetInt(id_get_exam_id)
#define id_new_patient                   930      // New patient command; val = 0;
#define id_new_exam                      238      // Start new exam for current
patient; val = 0;
#define id_quick_save                    107      // Image quick save; val = 0;
#define id_last_used_action_image_save    160      // Image save (save as) button's
command; val = 0;
#define id_quick_print                   108      // Image quick print; val = 0;
#define id_last_used_action_image_print    161      // Image printing button's
command; val = 0;
#define id_report_quick_print             135      // Report quick print; val = 0;
#define id_last_used_action_report        162      // Report button's command; val
= 0;
#define id_tool_measurement_1d            136      // Select "line" tool (primary
tool) of current mode; val = 0;
#define id_tool_measurement_2d            137      // Select "ellipse" tool
(secondary tool) of current mode; val = 0;
#define id_tool_annot_label               219      // Open Annotations ControlBar
and select label tool; val = 0;
#define id_tool_annot_arrow               220      // Open Annotations ControlBar
and select arrow tool; val = 0;
#define id_tool_annot_text                 237      // Open annotations ControlBar
and select text tool; val = 0;
#define id_tool_obj_delete                230      // "Delete selected" button; val
= 0;
#define id_measurements_delete_all         138      // "Delete All" button; val = 0;
#define id_thumbnails_show                635      // Show thumbnails; val = 0;
#define id_thumbnails_show_hide           638      // Show/hide thumbnails
depending on their current state; val = 0;
#define id_thumbnails_selection_shift      636      // Virtual marker for selecting
thumbnail for opening.
#define id_thumbnails_selection_open       637      // Open thumbnail that was
selected using virtual marker; val = 0;

```

```

#define id_presets                117    // Show Presets ControlBar; val
= 0;
#define id_preset_apply1          919    // Apply preset. As command
value is passed preset name. ParamSetString(id_preset_apply1, "Default")
#define id_presets_get_count      920    // GET the number of presets of
current probe. ParamGetInt(id_presets_get_count)
#define id_preset_get_name        921    // GET the name of preset of
currrent probe by passed preset index. ParamGetString2(id_preset_get_name, 0)
#define id_presets_get_count2     922    // GET the number of presets of
probe with passed code. ParamGetInt2(id_presets_get_count2, probe_code)
#define id_preset_get_name2       923    // GET the name of preset of
probe with passed code by passed preset index.
ParamGetString3(id_preset_get_name, probe_code, 0)
#define id_presets_load_from_path  933    // Load presets from passed
folder, for example, c:\temp\presets2. ParamSetString(id_presets_load_from_path,
"c:\\temp\\presets2")
#define id_cine_controls          246    // Show Cine ControlBar; val =
0;
#define id_measurements_controls  123    // Show Measurements ControlBar;
val = 0;
#define id_calculations_controls  124    // Show Calculations ControlBar;
val = 0;
#define id_navigate_calc_field    232    // When is opened Calculations
ControlBar, select previous/next field in it.
#define id_biopsy_controls        129    // Show Biopsy ControlBar; val =
0;
#define id_annotations_controls   130    // Show Annotations ControlBar;
val = 0;
#define id_body_marks_controls    131    // Show Body Marks ControlBar;
val = 0;
#define id_body_marks_get_count   233    // GET the number of body mark
images. ParamGetInt(id_body_marks_get_count)
#define id_body_mark_set_by_idx   234    // Set Body Mark by passed its
index (0, 1, 2, ...)
#define id_body_mark_remove       235    // Remove Body Mark; val = 0;
#define id_user_manual            144    // Open user manual; val = 0;
#define id_app_set_topmost        924    // Set application topmost
status: 1 - on, 0 - off
#define id_app_set_visible        925    // Set application visible
status: 1 - on, 0 - off
#define id_app_set_window_state   926    // Set main window state: 0 -
normal, 1 - minimized
#define id_app_exit               927    // Exit software; val = 0;
#define id_set_patient_name       928    // Set passed string as patient
name. ParamSetString(id_set_patient_name, "Text")
#define id_set_patient_id         929    // Set passed string as patient
id. ParamSetString(id_set_patient_id, "Text")
#define id_shutdown               186    // Initiate computer shut down;
val = 0;
#define id_patient2               931    // Show patient window; val = 0;
#define id_options2               932    // Show options window; val = 0;
// IMPORTANT. Cine end time is generated at the software graphical user
interface (GUI) after finishing of GUI invoked freeze command. This means that
it may differ from real time when hardware scanning was finished.
#define id_get_cine_end_date_time_str 690 // Get cine end (freeze)
absolute date and time as string. Format: "yyyy.MM.dd HH:mm:ss.ffffff".

```

ParamGetString(id\_get\_end\_date\_time\_str).