



## UP-D897 Black & White A6 Digital USB Printer

The Sony UP-D897 thermal printer is the ideal choice for digital ultrasound systems. This monochrome black and white printer incorporates a USB 2.0 interface for easy integration to digital ultrasound systems. The UP-D897 printer packs a powerful set of features in a small package. The chassis is extremely small and lightweight with a remarkable print speed of approximately 2.0 seconds.

- Ideal for digital ultrasound printing
- Smaller and lighter by 12% compared to the 895 series printer
- User-friendly controls on front panel
- Fast print speed — approximately 2.0 seconds in high speed mode
- USB 2.0 High Speed USB Interface easily connects to digital ultrasound system
- Uses the same media as the UP-895 series printers: UPP-110S/5, UPP-110HD, UPP-110HG



## Specifications

Printing Method	Direct Thermal Printing
Thermal Head	Thin film thermal 1280 dots
Resolution	325 dpi
Gradations	True 256 gray balance
Effective Print Pixels	4096 x 1280 (max.)
Interface	USB 2.0 High Speed Interface
Memory	4096 x 1280 x 8 bit
Print Speed	Approximately 2.0 seconds (High Speed Mode)
Dimensions	154mm x 88mm x 240mm (WHD) ( 6 1/8 x 3 1/2 x 9 1/2 inches - WHD)
Weight	Approximately 2.6kg (5lb. 12 oz)
Safety Standards	UL60950-1 Listed + CSA C22.2 No. 60950-1 + UL60601-1 Listed + CSA No. 60601-1
Supplied Accessories	AC Power Cord Operation Manual (CD ROM) Thermal Head Cleaning Sheet Driver Software UPP-110HG Media (1 Roll) USB Cable
Media (Not Included)	UPP-110HG High Gloss Paper UPP-110HD High Density Paper UPP-110S/5 Standard Paper
Driver Requirements	OS: Windows® 2000/XP



UP-D897 Rear Panel

## Print Media



UPP-110S/5  
Standard Media

UPP-110HD  
High Density Media

UPP-110HG  
Premium High Gloss Media

# SONY®

Sony Electronics Inc.  
1 Sony Drive  
Park Ridge, NJ 07656

[www.sony.com/medical](http://www.sony.com/medical)

©2007 Sony Electronics Inc. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features and specifications are subject to change without notice. All non-metric weights and measures are approximate. Sony is a trademark of Sony. Windows is a trademark of Microsoft. Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or other appropriately licensed medical professional. Caution: See product labeling for indications, contraindications, warnings, cautions and directions for use.

ME-UPD897