





SONDE CONVEX		Parametri				Applicazioni										
Foto	CODICE / Modello	Elementi Cristalli	Range Frequenze MHz	Raggio Curv. mm	Campo di Vista gradi mm	Abdominal	OB.-Gyn.	Vascular	Small Parts	MSK/Neuroi.	Pediatry	Cardiology	Transrectal	Transvaginal	Vascular Access	Interventional
 Convex HD	C5-2H60A-5	192 Elementi	2.0 - 5.0 MHz	65 mm	59°	■	■	■			■					

\* Our products are under constant development, so their design and specifications are indicative and may change at any time.

SONDE PHASED ARRAY		Parametri				Applicazioni										
Foto	CODICE / Modello	Elementi Cristalli	Range Frequenze MHz	Raggio Curv. mm	Campo di Vista Gradi mm	Abdomen	OB.-Gyn.	Vascular	Small Parts	MSK/Neuroi.	Pediatry	Cardiology	Transrectal	Transvaginal	Vascular Access	Interventional
 Phased Array Sector	P5-1S15-A6	64 Elementi	1.5 - 5.0 MHz	-	90°	■						■				

\* Our products are under constant development, so their design and specifications are indicative and may change at any time.



**TELEMED Medical Systems** Via  
Eugenio Villoresi, 24  
20143 Milano, Italy  
Phone: +39 02 36594100  
Mobile: +39 348 3190513  
Email: [info@telemedultrasound.com](mailto:info@telemedultrasound.com)  
Web: [www.telemedultrasound.com](http://www.telemedultrasound.com)

**TELEMED Ltd**  
Highway Business Center,  
Savanorių pr. 178A,  
LT-03154 Vilnius, Lithuania